

# NURA

## AI Early Detection Assistant

Detect earlier. Treat faster. Save lives.



gci solutions

Developed by GCT Solutions

# AGENDA

# NURA




- 1 The Problem & Market Opportunity
- 2 NURA Solution Overview
- 3 How NURA Works
- 4 Core Module 1: BEFAST (Stroke Early Detection)
- 5 Core Module 2: NIHSS (Stroke Severity Assessment)
- 6 Core Module 3: Autism Screening Mini-Game
- 7 Supporting Module: Health Tracking
- 8 Technology, Experience & Deployment


# THE COST OF LATE DETECTION

The challenge is not diagnosis-it's the lack of scalable, early, and standardized screening.

## STROKE: THE GOLDEN HOUR




## AUTISM: EARLY INTERVENTION WINDOW



# THE OPPORTUNITY

## Healthcare Is Shifting from Treatment to Early Detection

The demand for early screening is rising, and AI makes it possible - yet hospitals struggle to keep up




The opportunity is to deliver **fast, standardized, AI-assisted screening** - that scales with the needs of modern healthcare.

# Meet NURA

**Fast, accessible health screening at the point of care.**

An AI-powered self-screening assistant that provides quick, structured insights without staff dependency.



# The NURA Workflow

NURA enables fast, structured screening by combining user interaction with real-time AI analysis - across different care settings.

## Step 1: Access

User selects screening module (Mobile, Tablet, Kiosk).



## Step 2: Guided Interaction

Simple step-by-step instructions (movement, speech, or game). Takes 2-5 minutes.



## Step 3: AI Analysis

Real-time processing of visual speech, and behavioral signals



## Step 4: Result

Generation Structured outputs, risk indications, and supporting data.



## Step 5: Clinical

Review Results passed to clinicians for faster, informed decision-making.



# NURA - Main Module

From Early Screening to Clinical Decision - Powered by AI



## BEFAST (Stroke Early Detection)

A widely used method to highlight early warning signs of stroke (Balance, Eyes, Face, Arm, Speech, Time).

- **Users:** Public (patients, caregivers)
- **Devices:** Mobile, Kiosk
- **Context:** Hospital, public, at-home
- **Role:** Enables fast, self-guided screening to support early awareness and timely follow-up



## NIHSS (Stroke Severity Assessment)

A clinical scale used to evaluate stroke severity through structured neurological assessment.

- **Users:** Clinicians
- **Devices:** Mobile, Tablet
- **Context:** Hospital
- **Role:** Provides standardized scoring (0–42) to support treatment decisions after stroke is suspected



## AUTISM SCREENING

A brief screening to identify early behavioral and communication signs in children.

- **Users:** Children (with caregiver support)
- **Devices:** Mobile, Tablet
- **Context:** Home, school, controlled environments
- **Focus:** Attention, eye contact, and response patterns
- **Role:** Enables early behavioral signs to support timely




## Health Monitoring (Supporting Module)

Continuous tracking of key health indicators (BP, glucose, symptoms).

- **Users:** Patients / General users
- **Devices:** Mobile or integrated into hospital applications
- **Context:** Daily monitoring and clinical follow-up
- **Role:** Enables structured health tracking to support better clinical decisions over time

# CORE MODULE 1

# BEFAST




# Using AI to Support Screening for Early Warning Signs of Stroke with the BEFAST Framework

- **Built on the BEFAST Protocol:** A widely used clinical approach to identify possible stroke signs.
- **How AI Makes a Difference**
  - Faster: Quick assessment in minutes
  - More Consistent: Same standard applied every time
  - More Objective: Data-driven results to support decisions

## 6 test included:

- Facial Symmetry Check
- Arm Strength Test
- Leg Strength Test
- Balance Assessment
- Speech & Cognitive Evaluation
- Vision / Eye Tracking



## Core Module 1: BEFAST

# Facial Symmetry Check

### How the test works


- The patient looks directly at the screen
- The system guides the user step by step with voice instructions
- The patient performs simple actions:
  - Keep a neutral face (3 seconds)
  - Smile and hold (3 seconds)
  - Raise eyebrows and hold (3 seconds)

### What AI analyzes

- Compares left and right facial movement
- Detects asymmetry during expressions

### What the result shows

- Normal symmetry vs. potential facial drooping
- Indication of potential warning signs



## Core Module 1: BEFAST

# Arm Strength Test

### How the test works

- The patient raises both arms in front of the body
- Holds position for 10 seconds
- Follows guided instructions on screen

### What AI analyzes

- Tracks arm position over time
- Detects if one arm drifts downward

### What the result shows

- Stability vs. weakness
- Differences between left and right arms



## Core Module 1: BEFAST

# Leg Strength Test

### How the test works


- The patient lies or positions body ~1.5–2 meters from the screen
- Lifts each leg to ~30 degrees
- Holds for 5 seconds, then switches legs

### What AI analyzes

- Tracks leg height and stability
- Detects dropping or inability to hold position

### What the result shows

- Muscle strength and control
- Signs of weakness or imbalance



## Core Module 1: BEFAST

# Balance Assessment

### How the test works

- The patient stands following guided posture
- Holds balance position for a few seconds
- Repeats on both sides

### What AI analyzes

- Body posture and center of gravity
- Stability and body sway

### What the result shows

- Stable vs. unstable balance
- Signs of coordination issues

**Balance Analysis**

- Shoulder Offset: 2.21 mm ✓
- Hip Offset: 0.04 mm ✓
- Knee Support: 0.12 mm ✓
- Sway X: 0.00 mm ✓
- Sway Y: 0.12 mm ✓
- Leg Contact: 0.04 mm ✓
- Knee Flexion: 0.12 mm ✓

**Kiểm tra thăng bằng**

Điểm số hiện tại: **95**

Mức độ rủi ro: **Nguy cơ thấp**

**Chỉ số chi tiết**

Khung sườn lệch tư: 1.1 cm	Số lần đầu chỉnh thăng bằng: 0	Độ lệch trung bình: 1.1 cm
Thời gian kiểm tra: 20 s	Điểm hình ảnh: 0/4 (Trái: 0/4, Phải: 0/4)	Ngày mở đợt học: 0 %
Điểm bài kiểm tra: 98% - Xuất sắc - Giữ tư thế hoàn hảo	Điểm bài kiểm tra: 93% - Xuất sắc - Giữ tư thế hoàn hảo	Điểm bài kiểm tra: 3.0 %
Độ lệch tư đầu gối trái: 0.9 cm	Độ lệch tư đầu gối phải: 1.4 cm	

Chân trái

Giữ chân trái nâng lên. Giữ thăng bằng trên chân phải.  
Căng chân trái lên 90 độ với đầu gối gấp. Giữ chân phải thẳng và ổn định.

Bỏ qua Kết thúc

## Core Module 1: BEFAST

# Speech & Cognitive Test

### How the test works

- The patient follows guided tasks:
  - Describe a situation in an image
  - Name objects
  - Read sentences clearly
- The system records responses via microphone

### What AI analyzes

- Speech clarity
- Response accuracy
- Language and comprehension

### What the result shows

- Slurred or unclear speech
- Cognitive or language difficulties

The screenshot displays the BEFAST app interface. At the top, it shows the user's name 'Tùng Nguyễn' and the test title 'Kiểm tra lời nói'. A 'Speech Analysis' panel on the left lists metrics: Mouth Open (0.27), Mouth Vel (0.00), Mouth Base (0.00), Smile Alarm (0.00), and Lip Height Alarm (0.21). The main task area shows an illustration of a person painting a wall in a room with a dog and a doghouse. The instruction reads: 'Mô tả một tình huống bạn thấy trong hình ảnh. Nhấn Ghi âm khi sẵn sàng.' Below the illustration is a recording button and a '100%' completion indicator. A 'Kết quả kiểm tra lời nói' (Speech Test Results) window is overlaid, showing a score of 1 and a 'Không có bài ghi' (No recording) message. The results window includes a table of step scores and a 'Kết quả từng bước' (Step-by-step results) section.

Chi số chi tiết	Điểm
Bước 1: Đọc nội dung câu hỏi	1
Bước 2: Mô tả tình huống	Điểm: 1. Mô tả cơ bản nhưng thiếu một số chi tiết quan trọng và giảm độ trôi chảy.
Bước 3: Đọc tên đồ vật	Điểm: 1. Bạn đã mô tả đúng một trong các đồ vật nhưng không trôi chảy và thiếu thông tin về các đồ vật khác trong bức tranh.
Bước 4: Đọc câu lệnh	Điểm: 0. Lời nói rõ ràng, không có lỗi và dễ hiểu.
Bước 5: Đọc câu lệnh	Điểm: 0. Bệnh nhân đọc đúng câu mà không có lỗi, thể hiện khả năng giao tiếp tốt và hiểu biết rõ ràng.
Bước 6: Đọc câu lệnh	Điểm: 0. Lời nói rõ ràng, không có lỗi và dễ hiểu.
Thời gian hoàn thành	0 s

**Kết quả từng bước**

Bước	Điểm	Trạng thái
Bước 1: Đọc nội dung câu hỏi	1	Đạt
Bước 2: Mô tả tình huống	1	Đạt

## Core Module 1: BEFAST

# Vision/ Eye Tracking Test

### How the test works

- The patient follows a moving dot on the screen
  - Horizontal movement
  - Rectangular pattern
- Eye movement is continuously tracked

### What AI analyzes

- Gaze direction and tracking accuracy
- Focus and response consistency

### What the result shows

- Eye movement patterns
- Signs of visual or neurological issues

The image displays a patient undergoing an eye tracking test. The patient's eyes are tracked, with a red dot on the nose and green dots on the eyes. The background shows a computer screen with an 'Eye Analysis' panel and a sidebar menu. The 'Eye Analysis' panel lists various metrics with checkmarks, indicating successful completion. The sidebar menu includes options like 'Thông tin', 'Thời kỳ', 'Lịch sử', 'Kiểm tra', and 'Cài đặt'. The patient is sitting in front of a computer monitor displaying the test results.

**Kiểm tra chuyển động mắt**

Độ chính xác theo dõi tổng thể: **100 %**

Mức độ rủi ro: **Nguy cơ thấp**

Điểm NIHSS: **0 / 4** (Nguy cơ thấp)

KẾT QUẢ TUNG BƯỚC:

Nhìn sót nhất (Ngang)	0/1
Thị trường (Hình chữ nhật)	0/1
Đánh giá thị lực	0/2

Nguy cơ đột quỵ: **0.0 %**


Mình ảnh trước màn

Lưu | Làm lại kiểm tra | Đóng

## Core Module 1: BEFAST

# The Invisible Lens: Converting everyday patient behavior into clinical signals.

The system turns everyday patient actions into clear clinical information, helping clinicians make faster and more confident decisions.



# From Manual Screening to AI-Supported Workflow



## **BEFAST-Based Self-Screening (Entry Layer)**

Start screening on kiosk / mobile  
Perform guided BEFAST tasks with AI assistance  
AI analyzes facial, speech, and movement signals  
Automated potential warning signal



## **BEFAST (Self-Screening)**

At home: User searches for BEFAST guidance and performs self-checks  
At hospital: Doctor performs rapid BEFAST screening  
Based on visual observation and subjective judgment

## Impact on Early Stroke Detection



### For Healthcare Providers

- Earlier identification of potential stroke cases before clinical evaluation
- Reduced burden on initial screening and triage
- More consistent and standardized early screening



### For Patients

- Early awareness of stroke warning signs
- Faster action when symptoms appear
- Increased likelihood of timely medical attention



### System-Level Impact

- Scalable early screening beyond hospital settings
- Increased detection coverage across populations
- Reduced risk of missed early warning signs

# CORE MODULE 2

# NIHS

# S




## Core Module 2: NIHSS

# Structuring first-line stroke assessment with an AI Physician Assistant

Assisting clinicians with guided, step-by-step evaluations aligned with NIHSS protocols to determine stroke severity.


- ✔ **Structured clinical support:** Assists clinicians during patient evaluation with a guided approach
- ✔ **Post-screening assistance:** Supports assessment after initial screening signals
- ✔ **Standardized guidance:** Provides step-by-step evaluation aligned with NIHSS protocols
- ✔ **Consistent scoring support:** Enables structured scoring and reduces manual variation
- ✔ **Assist, not replace:** Supports clinical judgment without replacing clinician decisions



## Core Module 2: NIHSS

# Using AI to Support First-Line Stroke Assessment

## HOW IT WORK



## WHAT IS ASSESSED

Consciousness	Vision	Motor function
Language	Cognition	

0-42 NIHSS Scoring Scale	0: No stroke symptoms 1-4: Minor stroke	5-15: Moderate stroke	16-20: Moderate to severe stroke	21-42: Severe stroke
<b>1a. Level of consciousness</b> 0 = Alert, keenly responsive 1 = Not alert, but arousable by minor stimulation 2 = Not alert, requires repeated stimulation 3 = Unresponsive or responds only with reflex	<b>4. Facial palsy</b> 0 = Normal symmetric movements 1 = Minor paralysis 2 = Partial paralysis 3 = Complete paralysis of one or both sides	<b>7. Limb ataxia</b> 0 = Absent 1 = Present in one limb 2 = Present in two limbs	<b>8. Sensory</b> 0 = Normal; no sensory loss 1 = Mild-to-moderate sensory loss 2 = Severe to total sensory loss	<b>9. Best language</b> 0 = No aphasia; normal 1 = Mild to moderate aphasia 2 = Severe aphasia 3 = Mute, global aphasia
<b>1b. Level of consciousness questions</b> 0 = Answers two questions correctly 1 = Answers one question correctly 2 = Answers neither question correctly	<b>5. Motor arm</b> 0 = No drift 1 = Drift 2 = Some effort against gravity 3 = No effort against gravity; limb falls 4 = No movement	<b>10. Dysarthria</b> 0 = Normal 1 = Mild to moderate dysarthria 2 = Severe dysarthria	<b>11. Extinction and inattention</b> 0 = No abnormality 1 = Visual, tactile, auditory, spatial, or personal inattention 2 = Profound hemi-inattention or extinction	<b>Total score = 0-42</b>
<b>1e. Level of consciousness commands</b> 0 = Performs both tasks correctly 1 = Performs one task correctly 2 = Performs neither task correctly	<b>6. Motor leg</b> 6a. Left      6b. Right 0 = No drift 1 = Drift 2 = Some effort against gravity 3 = No effort against gravity 4 = No movement	<b>3. Visual</b> 0 = No visual loss 1 = Partial hemianopia 2 = Complete hemianopia 3 = Bilateral hemianopia		
<b>2. Best gaze</b> 0 = Normal 1 = Partial gaze palsy 2 = Forced deviation	<b>7. Limb ataxia</b> 0 = Absent 1 = Present in one limb 2 = Present in two limbs			

## Core Module 2: NIHSS

# From Manual Scoring to Standardized Stroke Severity Assessment

### 1. NIHSS (Today)



**OBSERVE & QUESTION:**  
Conducts NIHSS using checklist.



**NOTE ANSWERS:**  
Asks structured questions.



**CALCULATE SCORE:**  
Manually assigns score.



**Results vary by clinician and experience.**


### 2. AI-Assisted NIHSS (NURA)



**INITIATE ASSESSMENT:**  
Start guided NIHSS evaluation.



**AUTO INTERVIEW & ANALYSIS:**  
AI analyzes responses in real time.



### Clinical Impact

- ✓ Reduce variability in stroke severity assessment
- ✓ Support faster treatment decisions
- ✓ Improve consistency across departments
- ✓ Improve consistency across departments

**NURA delivers standardized, consistent stroke assessments with AI support.**

# Impact on Clinical Stroke Assessment

Faster and more consistent assessment supports timely treatment - improving overall stroke care outcomes.



## For Healthcare Providers

- Faster and more structured severity assessment
- Reduced manual effort in scoring and evaluation
- More consistent application of clinical standards



## For Patients

- Faster clinical evaluation after suspected stroke
- More timely and appropriate treatment decisions
- Reduced delays in care pathways



## System-Level Impact

- Standardized assessment across clinicians and settings
- Improved efficiency in stroke management workflows
- Better utilization of clinical resources


# CORE MODULE 3

# AUTISM DETECTION MINI-GAME



## What the Child Sees

An engaging mini-game with no formal instructions. The experience feels entirely like play.



## What the AI Sees

Continuous tracking of eye movement, attention/focus patterns, reaction timing, and behavioral consistency.

## Making autism screening engaging, scalable, and non-invasive

Game-based screening designed for non-clinical settings.  
It feels like play, but acts like a test.

**1** **Eye Movement & Gaze:**  
Tracking focus patterns.

**3** **Reaction Timing:**  
Measuring consistency and response speed.




**2** **Attention & Focus:**  
Analyzing engagement levels.

**4** **Behavioral Responses:**  
Recognizing patterns linked to autism indicators.

Generates an **Early Autism Risk Indication** and behavioral insights  
for specialist referral, without requiring a clinical setup.


## Core Module 3: Autism Detection

# Making Autism Screening Engaging, Scalable, and Non-Invasive



KEY BENEFITS	NUR	For Caregivers & Providers
<ul style="list-style-type: none"><li>✓ Non-invasive and child-friendly.</li><li>✓ Does not require clinical setup.</li></ul>	Enables earlier and more accessible autism screening without adding clinical complexity.	<ul style="list-style-type: none"><li>✓ Early signals for further evaluation.</li><li>✓ Better support for specialist referral decisions.</li></ul>

# SUPPORTING MODULE: HEALTH MONITORING



Eliminates memory bias for patients.



Provides clinicians with longitudinal data for faster, informed consultations.



## Supporting Module

# A Structured Health Journal for Better Clinical Review

### WHAT USERS LOG

- Daily Measurements: Blood Pressure, Blood Glucose
- Symptoms & Pain: Time, location, severity, notes
- Flexible Input: Voice recording or manual typing



### HOW NURA WORKS

- Capture & Timestamp: Records inputs with date and time
- AI Follow-up Questions: Asks simple questions to clarify details
- Structured Health Records: Organizes data consistently by type
- Trend Aggregation: Generates summaries and trends over time

### WHAT THIS ENABLES


- Accurate History Without Memory Bias: No need to rely on memory when visiting the doctor.
- Clear Trends Over Time: See patterns of pain, BP, and glucose across days and weeks.
- Ready-to-Share Report for Doctors: Export a monthly summary in one tap.
- Faster & More Informed Consultation: Doctors can assess your condition quickly and make better clinical decisions.

*From scattered inputs to structured health records - ready for clinical use.*



## Supporting Module

# From Daily Tracking to Clinical Insight



### For Patients

- Easier tracking without relying on memory
- Better visibility into health patterns

### For Clinicians

- Faster and more informed consultations
- Access to longitudinal patient data
- Improved decision-making

### Deployment

- Mobile-first experience
- Can be integrated into hospital systems

This module complements screening by enabling continuous health visibility.

**TECHNOLOGY**

**EXPERIENCE &  
DEPLOYMENT**



## User Experience

# Designed for absolute adoption, not complexity.


Simple for patients to use independently; valuable for clinicians making fast decisions.

## The Patient POV

- ✓ Takes only 2-5 minutes per check.
- ✓ Simple, self-guided step-by-step interaction.
- ✓ Zero staff required during screening.
- ✓ Available in multiple languages.



## The Doctor POV



Clear results available immediately.


Real-time dashboard with structured risk levels.

Easy to review and act on, with AI-supported insights

## AI & Technology Layer

# Built for accurate, real-world clinical use.

Advanced behavioral recognition paired with absolute data privacy.



## Accuracy & Validation

# Built for Accurate, Fast, and Real-World Use

Flexible infrastructure that integrates seamlessly across any care setting.



### The AI Layer

- Computer Vision (Face & Eye tracking)
- Speech Analysis
- Behavioral Pattern Recognition



### Privacy & Speed

On-Device Processing: Data is processed locally on the hardware, ensuring maximum patient privacy and reliability without cloud latency.



### Flexible Deployment

Deployable in hospitals, public screening locations, schools, or homes with no complex setup required.

Note: Detailed validation metrics available upon request / pilot deployment.

## Deployment Model

# Flexible deployment across the entire healthcare ecosystem.

Easy to deploy with zero complex setup required, bridging the gap between home and hospital.



### The Kiosk

Hospitals, Clinics, and  
Public Screening Locations.



### The Tablet

Dedicated screening stations  
and school settings.



### The Mobile Phone

Home and community  
tracking.

## Competitive Advantage

# What Makes NURA Different in Real-World Screening?

More than just a tool-it enables healthcare teams to act earlier, scale wider, and perform more consistently.

	Traditional Screening	NURA AI
Workflow	✓ Manual, staff-dependent.	✓ Self-service and automated.
Scope	✓ Single condition per assessment.	✓ Multi-condition (Stroke, Autism, Health tracking).
Speed	✓ Results delayed by clinical bottlenecks.	✓ Real-time results in 2-5 minutes.
Privacy	✗ Often cloud-based, data transferred externally.	✓ On-device processing for absolute data privacy.
Accessibility	✗ Strictly limited to specific medical facilities.	✓ Flexible across hospital, public, and home settings.

# Scaling early detection. Empowering clinicians. Transforming outcomes.

NURA makes vital early warning signs visible, standardized, and actionable.



System-level impact on overloaded hospital workflows.



Consistent, objective clinical data for every patient.



Capturing the intervention window before it closes.